

Tutela V4 - Release Notes

<u>Version</u>	<u>Description</u>	<u>Release Date</u>
01	Issue	15/06/2016
1.1	Issue	29/09/2016

Contents

1.0 Overview	2
2.0 Additional Features	2
4.0 Removed/Retired Features	2
5.0 Bug Fixes	2
6.0 System Requirements	3

1.0 Overview

1.1 This document outlines the additional features, Improvements and bug fixes to the 4th generation website.

1.2 This version as with all previous versions provides remote access to the automatic recorded data and alarm exceptions of the sensing probes installed in the customer's facility. Customers can gain access to their data via the secure Tutela V4 web base application hosted by Tutela medical. The data records are published in real time via these secure web pages. Remote site hosting, maintenance and validation of electronic data records are regulated through Tutela Medical **UK ISO 9001 Business Management System**.

2.0 Additional Features

2.1 Access to the accelerated ageing application and library for selected users.

2.2 Help Functionality added to the User Management.

2.3 Hit counter displayed to users with relevant permission.

3.0 Improvements/Enhancements

3.1 Inches water gauge displaying to 4 decimal places. (CR-012)

3.2 The automatic graph scaling function is now selectable.

3.3 Pin codes can be entered by keyboard in addition to the drop down.

3.4 The check readings functionality timeout has been reduced enabling the user to advance in a timelier fashion if all sensor readings are not retrieved.

3.5 Improved response time for sensor value loading.

4.0 Removed/Retired Features

4.1 Automatic auto scaling on sensor graphs.

5.0 Bug Fixes

5.1 User management division selection showing retired customers.

- 5.2 No pin code requested for uploading or creating documents to pdf library.
- 5.3 Alarm contact permission when not allowing access.
- 5.4 Users without permission accessing setting page though bottom side menu.
- 5.5 Old technology not showing time delay on setting page.

6.0 System Requirements

6.1 Operating Systems and browsers compatibility is as follows:

Yes
No

	Windows 10	Windows 7 64-Bit	Windows 7	Windows XP
Edge		N/A	N/A	N/A
IE 11			N/A	N/A
IE 10	N/A		N/A	N/A
IE 9 x 64	N/A			N/A
IE 9	N/A			N/A
IE 8	N/A			
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48		N/A		
Chrome 47		N/A		
Chrome 46		N/A		
Chrome 45		N/A		
Chrome 44		N/A		
Chrome 43		N/A		
Chrome 42	N/A	N/A		
Chrome 41	N/A	N/A		
FireFox 47				
FireFox 46				
FireFox 45				
Firefox 44	N/A	N/A		
Firefox 43	N/A	N/A		
Firefox 42	N/A	N/A		
Firefox 41	N/A	N/A	N/A	N/A
Firefox 40	N/A	N/A	N/A	N/A
Firefox 39	N/A	N/A	N/A	N/A
Firefox 38	N/A	N/A	N/A	N/A

Opera 37			N/A	N/A
Opera 36				
Opera 35				
Safari 9	N/A	N/A	N/A	N/A
Safari 8	N/A	N/A	N/A	N/A
Safari 7	N/A	N/A	N/A	N/A
Safari 6.2	N/A	N/A	N/A	N/A

	MAC 10.11	MAC 10.10	MAC 10.9	MAC 10.8
Edge	N/A	N/A	N/A	N/A
IE 11	N/A	N/A	N/A	N/A
IE 10	N/A	N/A	N/A	N/A
IE 9 x 64	N/A	N/A	N/A	N/A
IE 9	N/A	N/A	N/A	N/A
IE 8	N/A	N/A	N/A	N/A
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48				
Chrome 47				
Chrome 46				
Chrome 45				
Chrome 44				
Chrome 43				
Chrome 42				
Chrome 41				
FireFox 47				
FireFox 46				
FireFox 45				
Firefox 44				
Firefox 43				
Firefox 42				
Firefox 41	N/A	N/A	N/A	N/A
Firefox 40	N/A	N/A	N/A	N/A
Firefox 39	N/A	N/A	N/A	N/A
Firefox 38	N/A	N/A	N/A	N/A
Opera 37				N/A
Opera 36				
Opera 35				
Safari 9			N/A	N/A
Safari 8	N/A		N/A	N/A
Safari 7	N/A	N/A		N/A
Safari 6.2	N/A	N/A	N/A	

