

Tutela V4.1.5 release notes



Version	Description	Release Date
0.1	Issue	15/06/2016
1.1	Issue	29/09/2016
1.2	Issue	03/02/2017
1.3	Issue	04/07/2017
1.4	Issue	03/11/2017
1.5	Issue	21/05/2018

# contents

1.0 Overview	2	
2.0 Additional Features	2	
3.0 Improvements/Enhancements	3	
4.0 Removed/Retired Features	3	
5.0 Bug Fixes	3	
6.0 System Requirements	4	







### 1.0 Overview

- **1.1** This document outlines the additional features, Improvements, and bug fixes to the 4th generation interface, version 4.1.4.
- 1.2 This version as with all previous versions provides remote access to the automatic recorded data and alarm exceptions of the sensing probes installed in the customer's facility. Customers can gain access to their data via the secure Tutela V4 web base application hosted by Tutela medical. The data records are published in real time via these secure web pages. Remote site hosting, maintenance and validation of electronic data records are regulated through Tutela Medical UK ISO 9001 Business Management System.

### 2.0 Additional Features

- **2.1** Improved login security Bcrypt encryption implemented for tighter password security.
- 2.2 Training Videos Training videos are now available to all users to aid user experience (CR 190).









## 3.0 Improvements/Enhancements

- **3.1** WebGui Navigation Any Websites supporting WebGui access will navigate back to Tutela v4.
- **3.2** Login The interface now allows for special characters in the user name.
- **3.3** Library Reorganization of the year folders so all old documents created before the new folder system was implemented are contained within a folder named 2016-2017.
- **3.4** Incident Reports Special characters now allowed in the input fields of the incident reports.
- **3.5** Library Pending file directory now allows web administrators to select multiple files to unload.
- **3.6** Graphing Improvements to accuracy of the incidents showing on the graph.
- **3.7** Mkt Calculator Precision has increased to 4 decimal places.

#### 4.0 Removed/Retired Features

**4.1** None.

## 5.0 Bug Fixes

- **5.1** Reporting Pdf not allowing secondary signoff when using internet explorer version 11.0.28 (windows 7). A download button has now been implemented to allow these users to perform secondary signoff's.
- **5.2** Mkt Calculator This function has been adjustment to perform its calculation gas constant of DELTA\_H\_OVER\_R=1000.
- **5.3** Sensor List Page Sorting the sensors stopped the sensor values from loading.







РС	Windows 10	Windows 7 64-Bit	Windows 7	Windows XP
Edge				
IE 11				
IE 10				
IE 9 x 64				
IE 9				
IE 8				
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48				
Chrome 47				
Chrome 46				
Chrome 45				
Chrome 44				
Chrome 43				
Chrome 42				
Chrome 41				
FireFox 47				
FireFox 46				
FireFox 45				
FireFox 44				
FireFox 43				
FireFox 42				
FireFox 41				
FireFox 40				





yes	no n/a			
	Windows 10	Windows 7 64-Bit	Windows 7	Windows XP
Firefox 39				
Firefox 38				
Opera 37				
Opera 36				
Opera 35				
Safari 9				
Safari 8				
Safari 7				
Safari 6.2				





yes no n/a

MAC	MAC 10.11	MAC 10.10	MAC 10.9	MAC 10.8
Edge				
IE 11				
IE 10				
IE 9 x 64				
IE 9				
IE 8				
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48				
Chrome 47				
Chrome 46				
Chrome 45				
Chrome 44				
Chrome 43				
Chrome 42				
Chrome 41				
FireFox 47				
FireFox 46				
FireFox 45				
FireFox 44				
FireFox 43				
FireFox 42				
FireFox 41				
FireFox 40				









