



WIRELESS TEMPERATURE MONITORING SYSTEM GENESIS 3

Tutela V4.1.6 - Release Notes

Tutela V4.1.6 - Release Notes

<u>Version</u>	<u>Description</u>	<u>Release Date</u>
01	Issue	15/06/2016
1.1	Issue	29/09/2016
1.2	Issue	03/02/2017
1.3	Issue	04/07/2017
1.4	Issue	03/11/2017
1.5	Issue	01/03/2018
1.6	Issue	10/10/2018

Contents


1.0 Overview	4
2.0 Additional Features	4
4.0 Removed/Retired Features	5
5.0 Bug Fixes	5
6.0 System Requirements	6

1.0 Overview

- 1.1 This document outlines the additional features, Improvements, and bug fixes to the 4th generation website, version 4.1.5.
- 1.2 This version as with all previous versions provides remote access to the automatic recorded data and alarm exceptions of the sensing probes installed in the customer’s facility. Customers can gain access to their data via the secure Tutela V4 website application hosted by Tutela medical. The data records are published in real time via these secure web pages. Remote site hosting, maintenance and validation of electronic data records are regulated through Tutela Medical **UK ISO 9001 Business Management System**.

2.0 Additional Features

- 2.1 Alarms Report Page – A monthly alarms report will be generated and made available to be viewed on the website. A user with the correct permissions will be able to download or save this report to their customer library. These reports will summarise the sites incident activity for the period. The report details alarms received, alarms actioned customer and Tutela response times and provides a breakdown of the activity. (CR 175)
- 2.2 Site List Page - PQ Web site access, a function to give Tutela support the ability to grant a 14-day temporary window which provides new customers a chance to train on the new system and produce their PQ’s (CR 200).



- 2.3 Isolation History – simplified table layout, making it easier to read with an addition of an early termination column, if the isolation is removed before is a natural end. (CR 256)
- 2.4 User password reminder service – A service to remind the user their password is about to expire. (CR 201)

3.0 Improvements/Enhancements

- 3.1 Compare Graphs –An increase to the size of ledger text. (CR 258)
- 3.2 Contacts - USA 001 prefix. When adding new USA contracts phone numbers to the website a 001 prefix will now be standard. (CR 177).
- 3.3 Incidents – The Tutela operator column has been replaced by the cause identified extracted from the incident report, giving better identification of incidents. (CR 255)
- 3.4 Library – Check reading directories will be sorted in descending date order.
- 3.5 Login – forgotten password will force user to reset after email is sent out.
- 3.6 Pending Files – The folder has been added at a site level, allowing more flexibility for customer admins so they can copy the support documentation to the relevant departments.

4.0 Removed/Retired Features

5.0 Bug Fixes

- 5.1 Compare Graphs Report – Null values caused creation of report to fail.
- 5.2 Compare Graphs Report – gaps in the graph data are now reflected.
- 5.3 Settings – Alarm parameters for inH2O can be set with a precision of 4 decimal places.

6.0 System Requirements

6.1 Operating Systems and browsers compatibility is as follows:

Yes
No

	Windows 10	Windows 7 64-Bit	Windows 7	Windows XP
Edge		N/A	N/A	N/A
IE 11			N/A	N/A
IE 10	N/A		N/A	N/A
IE 9 x 64	N/A			N/A
IE 9	N/A			N/A
IE 8	N/A			
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48		N/A		
Chrome 47		N/A		
Chrome 46		N/A		
Chrome 45		N/A		
Chrome 44		N/A		
Chrome 43		N/A		
Chrome 42	N/A	N/A		
Chrome 41	N/A	N/A		
FireFox 47				
FireFox 46				
FireFox 45				
Firefox 44	N/A	N/A		
Firefox 43	N/A	N/A		
Firefox 42	N/A	N/A		
Firefox 41	N/A	N/A	N/A	N/A
Firefox 40	N/A	N/A	N/A	N/A
Firefox 39	N/A	N/A	N/A	N/A
Firefox 38	N/A	N/A	N/A	N/A
Opera 37			N/A	N/A
Opera 36				
Opera 35				
Safari 9	N/A	N/A	N/A	N/A
Safari 8	N/A	N/A	N/A	N/A
Safari 7	N/A	N/A	N/A	N/A
Safari 6.2	N/A	N/A	N/A	N/A

	MAC 10.11	MAC 10.10	MAC 10.9	MAC 10.8
Edge	N/A	N/A	N/A	N/A
IE 11	N/A	N/A	N/A	N/A
IE 10	N/A	N/A	N/A	N/A
IE 9 x 64	N/A	N/A	N/A	N/A
IE 9	N/A	N/A	N/A	N/A
IE 8	N/A	N/A	N/A	N/A
Chrome 51				
Chrome 50				
Chrome 49				
Chrome 48				
Chrome 47				
Chrome 46				
Chrome 45				
Chrome 44				
Chrome 43				
Chrome 42				
Chrome 41				
FireFox 47				
FireFox 46				
FireFox 45				
Firefox 44				
Firefox 43				
Firefox 42				
Firefox 41	N/A	N/A	N/A	N/A
Firefox 40	N/A	N/A	N/A	N/A
Firefox 39	N/A	N/A	N/A	N/A
Firefox 38	N/A	N/A	N/A	N/A
Opera 37				N/A
Opera 36				
Opera 35				
Safari 9			N/A	N/A
Safari 8	N/A		N/A	N/A
Safari 7	N/A	N/A		N/A
Safari 6.2	N/A	N/A	N/A	